




**Product Specification**  
**Absorbent liner, suitable for mailing container 84  
x 30 mm, (LxW): 75 x 65 mm**

**Page 1**



**Product description**

Order number 78.684  
Product description Absorbent liner, suitable for mailing container 84 x 30 mm,  
(LxW): 75 x 65 mm, 100 piece(s)/case

**Product characteristics**

Application Mailing of samples, Multi-Safe disposel boxes  
Suitable for mailing container 84 x 30 mm

**Size**

Width of product 65 mm  
Length of product 75 mm

**Material & colours**

Product material cellulose

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Absorbent liner, suitable for mailing container 84**  
**x 30 mm, (LxW): 75 x 65 mm**

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      2000  
Type of smallest                              case  
subpackaging  
Piece(s) / inner box                      100  
Piece (s) / outer case                      2000  
Piece(s) / pallet                              120000  
Depth of case                              410 mm  
Width of case                                260 mm  
Height of case                              305 mm  
Case volume                                0.0325 cbm  
Weight of product                        0.0007 kg  
Weight of case                              1.9 kg  
EAN of inner box                        4038917083362  
EAN case                                    4038917002295

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

