




Product Specification
**Absorbent liner, suitable for mailing container 84
x 30 mm, (LxW): 75 x 65 mm**

Page 1



Product description

Order number	78.684
Product description	Absorbent liner, suitable for mailing container 84 x 30 mm, (LxW): 75 x 65 mm, 100 piece(s)/case

Product characteristics

Application	Mailing of samples, Multi-Safe disposel boxes
Suitable for	mailing container 84 x 30 mm

Size

Width of product	65 mm
Length of product	75 mm

Material & colours

Product material	cellulose
------------------	-----------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Absorbent liner, suitable for mailing container 84
x 30 mm, (LxW): 75 x 65 mm

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 2000
Type of smallest case
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 2000
Piece(s) / pallet 120000
Depth of case 410 mm
Width of case 260 mm
Height of case 305 mm
Case volume 0.0325 cbm
Weight of product 0.0007 kg
Weight of case 1.9 kg
EAN of inner box 4038917083362
EAN case 4038917002295

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

