



Product description

| | |
|---------------------|---|
| Order number | 14.1513.002 |
| Product description | Colour-coded cap, for S-Monovette® , yellow, 100 piece(s)/bag |

Product characteristics

| | |
|--------------|--------------|
| Suitable for | S-Monovette® |
|--------------|--------------|

Size

| | |
|-------------------|-------|
| Length of product | 13 mm |
|-------------------|-------|

Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | High Density Polyethylene (HD-PE) |
| Colour of product | yellow |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 324000 |
| Depth of case | 268 mm |
| Width of case | 191 mm |
| Height of case | 100 mm |
| Case volume | 0.0051 cbm |
| Weight of product | 0.0003 kg |
| Weight of case | 0.47 kg |
| EAN of inner box | 4038917251556 |
| EAN case | 4038917251549 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

