



### **Product description**

Order number	14.1513
Product description	Colour-coded cap, for S-Monovette® , pink, 100 piece(s)/bag

### **Product characteristics**

Suitable for	S-Monovette®
Connection	S-Monovette®

### **Size**

Length of product	13 mm
-------------------	-------

### **Material & colours**

Product material	High Density Polyethylene (HD-PE)
Colour of product	pink

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	324000
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.0003 kg
Weight of case	0.47 kg
EAN of inner box	4038917086790
EAN case	4038917073943

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

