



### **Product description**

Order number	53.1030
Product description	Plunger, translucent, for Seraplas® V filter, 200 piece(s)/bag

### **Product characteristics**

Suitable for	Seraplas® V filter
--------------	--------------------

### **Size**

Length of product	89 mm
-------------------	-------

### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	translucent

### **Purity & certification**

Product category	no medical device   IVD
Safety data sheet	no
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	200
Piece (s) / outer case	2000
Piece(s) / pallet	96000
Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0013 kg
Weight of case	3.12 kg
EAN of inner box	4038917085076
EAN case	4038917003773

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

