



**Product Specification**  
**Seraplas® valve filter, green, for separation of  
serum/plasma from the blood cells after  
centrifugation**

**Page 1**



**Product description**

Order number	53.420
Product description	Seraplas® valve filter, V16, green, for separation of serum/plasma from the blood cells after centrifugation, 100 piece(s)/case

**Product characteristics**

Application	venous blood collection
Suitable for	separation of serum/plasma from the blood cells after centrifugation

**Size**

Length of product	92 mm
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Seraplas® valve filter, green, for separation of serum/plasma from the blood cells after centrifugation**

**Material & colours**

Product material	Polyethylene (PE)
Colour of product	green

**Purity & certification**

Product category	no medical device   IVD
Safety data sheet	no
Batched	no

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	54000
Depth of box	156 mm
Width of box	99 mm
Height of box	130 mm
Depth of case	513 mm
Width of case	328 mm
Height of case	150 mm
Case volume	0.0252 cbm
Weight of product	0.0023 kg
Weight of case	2.98 kg
EAN of inner box	4038917081498
EAN case	4038917000581

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

