



### **Product description**

Order number	55.1571.051
Product description	Adapter tube, (LxØ): 65 x 13 mm, material: PP, light blue, 250 piece(s)/bag

### **Product characteristics**

Suitable for	S-Monovette® with Ø 8 mm, Microvette® and micro tubes
Base shape	round base

### **Size**

Diameter	13 mm
Length of product	65 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	light blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	1000
Piece(s) / pallet	72000
Depth of case	333 mm
Width of case	303 mm
Height of case	190 mm
Case volume	0.0192 cbm
Weight of product	0.0019 kg
Weight of case	2.23 kg
EAN of inner box	4038917283380
EAN case	4038917283373

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

