



Product description

Order number	55.1571
Product description	Adapter tube, (LxØ): 65 x 13 mm, material: PP, transparent, 250 piece(s)/bag

Product characteristics

Closure type	screw cap
Suitable for	S-Monovette® with Ø 8 mm, Microvette® and micro tubes
Base shape	round base

Size

Diameter	13 mm
Length of product	65 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1000
Type of smallest bag
subpackaging
Piece(s) / inner box 250
Piece (s) / outer case 1000
Piece(s) / pallet 72000
Depth of case 333 mm
Width of case 303 mm
Height of case 190 mm
Case volume 0.0192 cbm
Weight of product 0.0019 kg
Weight of case 2.23 kg
EAN of inner box 4038917242523
EAN case 4038917242516

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

