



Product description

Order number	55.525.002
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PP, round base, brown, push cap, suitable for Microvette® and micro tube with screw cap, 250 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Suitable for	Microvette® and micro tube with screw cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g

Size

Diameter	13 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tube, 5 ml, (LxØ): 75 x 13 mm, PP

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	brown

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Depth of case	398 mm
Width of case	398 mm
Height of case	265 mm
Case volume	0.042 cbm
Weight of product	0.0025 kg
Weight of case	5.57 kg
EAN of inner box	4038917109147
EAN case	4038917012355

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-02

This document was prepared by EDP support and is valid without signature.



SARSTEDT