



### **Product description**

Order number	57.477.500
Product description	Tube, working volume: 4.5 ml, (LxØ): 75 x 12 mm, material: PP, conical base, transparent, push cap, suitable for Pharmacia/LKB Gamma Counter and Gamma Counter Cobra from Packard, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Suitable for	Pharmacia/LKB Gamma Counter and Gamma Counter Cobra from Packard
Base shape	conical base
Volume of work	4,5 ml
Centrifugation max (RCF)	5000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Diameter	12 mm
Length of product	75 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	4000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	4000
Piece(s) / pallet	96000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0025 kg
Weight of case	10.59 kg
EAN of inner box	4038917107778
EAN case	4038917011945

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

