



Product description

| | |
|---------------------|---|
| Order number | 65.163 |
| Product description | Screw cap, material: HD-PE, natural, for tubes 75 x 13 mm, 1,000 piece(s)/bag |

Product characteristics

| | |
|--------------|------------------|
| Suitable for | tubes 75 x 13 mm |
|--------------|------------------|

Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | High Density Polyethylene (HD-PE) |
| Colour of product | natural |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 5000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 5000 |
| Piece(s) / pallet | 210000 |
| Depth of case | 398 mm |
| Width of case | 398 mm |
| Height of case | 250 mm |
| Case volume | 0.0396 cbm |
| Weight of product | 0.0014 kg |
| Weight of case | 7.45 kg |
| EAN of inner box | 4038917272223 |
| EAN case | 4038917272216 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

