



Product description

Order number	65.803
Product description	Push cap, natural, suitable for tubes Ø 15.7 mm & cuvette 67.743, 1,000 piece(s)/bag

Product characteristics

Cap	standard cap
Closure type	push cap
Suitable for	tubes Ø 15.7 mm & cuvette 67.743

Size

Diameter	16.35 mm
Height of product	20.45 mm

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	210000
Depth of case	398 mm
Width of case	398 mm
Height of case	250 mm
Case volume	0.0396 cbm
Weight of product	0.0013 kg
Weight of case	7.1 kg
EAN of inner box	4038917084192
EAN case	4038917003056

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

