



### **Product description**

Order number	65.806
Product description	Push cap, natural, suitable for tubes Ø 13 mm, 1,000 piece(s)/bag

### **Product characteristics**

Cap	standard cap
Closure type	push cap
Suitable for	tubes Ø 13 mm

### **Size**

Diameter	13.3 mm
Height of product	15.2 mm

### **Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      10000  
Type of smallest subpackaging                      bag  
Piece(s) / inner box                      1000  
Piece (s) / outer case                      10000  
Piece(s) / pallet                      420000  
Depth of case                      398 mm  
Width of case                      398 mm  
Height of case                      250 mm  
Case volume                      0.0396 cbm  
Weight of product                      0.0007 kg  
Weight of case                      7.36 kg  
EAN of inner box                      4038917084222  
EAN case                      4038917003087

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

