



Product description

Order number	95.1008
Product description	One-patient tourniquet, material: polymer, latex-free, 25 piece(s)/bag

Product characteristics

Application	venous blood collection
-------------	-------------------------

Material & colours

Product material	polymer
Colour of product	yellow

Purity & certification

Product category	Medical device
CE certified	CE supplier
Latex-free	yes
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	25
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	500
Piece(s) / pallet	8000
Depth of case	390 mm
Width of case	230 mm
Height of case	560 mm
Case volume	0.0502 cbm
Weight of product	0.0087 kg
Weight of case	5.2 kg
EAN of inner box	5060728810201
EAN case	5060728810218

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

