



Product description

| | |
|---------------------|--|
| Order number | 85.1160 |
| Product description | S-Monovette® needle, needle (ØxL): 20G x 1 1/2", coding colour: yellow, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 100 piece(s)/case |

Product characteristics

| | |
|----------------|--------------------|
| Coding colour | yellow |
| Safety product | no |
| Connection | membrane screw cap |

Size

| | |
|-------------------------------|---------------|
| Diameter of needle (external) | 20 G / 0,9 mm |
| Length of needle | 38 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**S-Monovette® needle, 20G x 1 1/2", yellow, 1
piece(s)/blister**

Page 2

Purity & certification

| | |
|------------------|---------------------------------------|
| Product category | Medical device |
| CE certified | CE 0124 (DEKRA) |
| Purity standard | sterile, non-pyrogenic/endotoxin-free |
| Sterilisation | Gamma irradiation |
| Pyrogen-free | yes |
| Batched | yes |

Packaging

| | |
|----------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 78000 |
| Depth of box | 268 mm |
| Width of box | 111 mm |
| Height of box | 71 mm |
| Depth of case | 373 mm |
| Width of case | 233 mm |
| Height of case | 290 mm |
| Case volume | 0.0252 cbm |
| Weight of product | 0.0018 kg |
| Weight of case | 2.96 kg |
| EAN of inner box | 4038917085519 |
| EAN case | 4038917004169 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-06

This document was prepared by EDP
support and is valid without signature.



SARSTEDT