



Product Specification

S-Monovette® safety needle, 21G x 1 1/2", green,
1 piece(s)/blister

Page 1



Product description

Order number	85.1162.200
Product description	S-Monovette® safety needle, needle (ØxL): 21G x 1 1/2", coding colour: green, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 50 piece(s)/case
PZN number	6796201

Product characteristics

Coding colour	green
Safety product	yes
Connection	membrane screw cap

Size

Diameter of needle (external)	21 G / 0,8 mm
Length of needle	38 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-08

This document was prepared by EDP
support and is valid without signature.



SARSTEDT



Product Specification

**S-Monovette® safety needle, 21G x 1 1/2", green,
1 piece(s)/blister**

Page 2

Purity & certification

Product category	Medical device
CE certified	CE 0124 (DEKRA)
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	blister
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	39000
Depth of box	268 mm
Width of box	111 mm
Height of box	71 mm
Depth of case	373 mm
Width of case	233 mm
Height of case	290 mm
Case volume	0.0252 cbm
Weight of product	0.0025 kg
Weight of case	2.45 kg
EAN of inner box	4038917126045
EAN case	4038917361927

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-08

This document was prepared by EDP support and is valid without signature.



SARSTEDT