



### **Product description**

Order number	85.1372.200
Product description	S-Monovette® safety needle, needle (ØxL): 20G x 1", coding colour: yellow, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 50 piece(s)/case

### **Product characteristics**

Coding colour	yellow
Safety product	yes
Connection	membrane screw cap

### **Size**

Diameter of needle (external)	20 G / 0,9 mm
Length of needle	25.4 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	Medical device
CE certified	CE 0124 (DEKRA)
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	blister
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	39000
Depth of box	268 mm
Width of box	111 mm
Height of box	71 mm
Depth of case	373 mm
Width of case	233 mm
Height of case	290 mm
Case volume	0.0252 cbm
Weight of product	0.0026 kg
Weight of case	2.47 kg
EAN of inner box	4038917127943
EAN case	4038917361941

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

