




**Product Specification**  
**Safety-Multifly® needle, 23G x 3/4", blue, tube**  
**length: 200 mm, 1 piece(s)/blister**

**Page 1**



**Product description**

Order number	85.1640.935
Product description	Safety-Multifly® needle, needle (ØxL): 23G x 3/4", coding colour: blue, tube length: 200 mm, with needle holder assembled, DEHP-free, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 25 piece(s)/case
PZN number	6797488

**Product characteristics**

Design	needle holder assembled
Coding colour	blue
Safety product	yes
Connection	membrane screw cap

**Size**

Diameter of needle (external)	23 G
Diameter of tube	2.3 mm
Length of needle	19 mm
Length of tube	200 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Safety-Multifly® needle, 23G x 3/4", blue, tube**  
**length: 200 mm, 1 piece(s)/blister**

**Page 2**

**Purity & certification**

Product category	Medical device
CE certified	CE 0124 (DEKRA)
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Plasticiser	DEHP-free
Batched	yes

**Packaging**

Minimum order qty.	10000
Type of smallest subpackaging	blister
Piece(s) / inner box	25
Piece (s) / outer case	100
Piece(s) / pallet	7200
Depth of case	533 mm
Width of case	263 mm
Height of case	140 mm
Case volume	0.0196 cbm
Weight of product	0.0052 kg
Weight of case	1.38 kg
EAN of inner box	4038917423410
EAN case	4038917423403

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-05

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**