



Product description

| | |
|---------------------|--|
| Order number | 85.925 |
| Product description | Micro needle, for collecting very small amounts of blood, needle (ØxL): 25G x 3/4", coding colour: orange, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 50 piece(s)/case |

Product characteristics

| | |
|----------------|--|
| Application | for collecting very small amounts of blood |
| Coding colour | orange |
| Safety product | no |

Size

| | |
|-------------------------------|---------------|
| Diameter of needle (external) | 25 G / 0,5 mm |
| Length of needle | 19 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

| | |
|------------------|---------------------------------------|
| Product category | Medical device |
| CE certified | CE 0124 (DEKRA) |
| Purity standard | sterile, non-pyrogenic/endotoxin-free |
| Sterilisation | Gamma irradiation |
| Pyrogen-free | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 50 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 55000 |
| Depth of box | 139 mm |
| Width of box | 108 mm |
| Height of box | 66 mm |
| Depth of case | 338 mm |
| Width of case | 228 mm |
| Height of case | 160 mm |
| Case volume | 0.0123 cbm |
| Weight of product | 0.0013 kg |
| Weight of case | 1.36 kg |
| EAN of inner box | 4038917281843 |
| EAN case | 4038917281836 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

