



Product Specification
Sample tube, Serum CAT, 4.5 ml, cap white,
(LxØ): 75 x 13 mm, with print

Page 1



Product description

Order number	32.329
Product description	Sample tube, nominal volume: 4.5 ml, cap white, preparation: Clotting activator, push cap, (LxØ): 75 x 13 mm, round base, with print, 100 piece(s)/bag

Product characteristics

Type of preparation	Clotting activator
Colour code	EU
Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	push cap
Base shape	round base
Recommended centrifugation	2000 x g
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Sample tube, Serum CAT, 4.5 ml, cap white,
(LxØ): 75 x 13 mm, with print

Page 2

Size

Sample volume	4.5 ml
Diameter	13 mm
Length of product	75 mm
Length excluding cap	75 mm

Material & colours

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	white
Tube material	Polypropylene (PP)

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	64000
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0.0255 cbm
Weight of product	0.0038 kg
Weight of case	4.23 kg
EAN of inner box	4038917081160
EAN case	4038917000307

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT