



Product description

Order number	41.1504.015
Product description	Micro sample tube EDTA K3E, venous blood collection, preparation: K3 EDTA, nominal volume: 1.3 ml, Ø: 10.8 mm, push cap, cap: red, colour code EU, with paper label, label/print: red, skirted conical base, 100 piece(s)/bag

Product characteristics

Type of collection	venous
Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	red
Graduation	yes
Cap	flat cap
Closure type	push cap
Application	venous blood collection
Base shape	skirted conical base
Volume of work	1,3 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Sample volume	1.3 ml
Diameter	10.8 mm
Length of product	44 mm
Length excluding cap	40 mm

Material & colours

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	red
Tube material	Polypropylene (PP)

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	96000
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0.0135 cbm
Weight of product	0.0016 kg
Weight of case	1.86 kg
EAN of inner box	4038917310222
EAN case	4038917310215

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

