



Product Specification
**Micro sample tube Fluoride/heparin FH, 1.3 ml,
push cap, EU**

Page 1



Product description

Order number	41.1505.005
Product description	Micro sample tube Fluoride/heparin FH, venous blood collection, preparation: Fluoride + heparin, nominal volume: 1.3 ml, (LxØ) with cap: 44 x 10.8 mm, push cap, cap: yellow, colour code EU, with plastic label, label/print: yellow, skirted conical base, 100 piece(s)/bag

Product characteristics

Type of collection	venous
Type of preparation	Fluoride + heparin
Preparation concentration	1.0 mg/ml blood + 16 IU/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	yellow
Graduation	yes
Cap	flat cap
Closure type	push cap
Application	venous blood collection
Base shape	skirted conical base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Micro sample tube Fluoride/heparin FH, 1.3 ml,
push cap, EU**

Page 2

Volume of work 1,3 ml
Recommended centrifugation 2000 x g
Centrifugation time 5 Min

Size

Sample volume 1.3 ml
Diameter 10.8 mm
Length of product 44 mm
Length including cap 44 mm
Length excluding cap 40 mm

Material & colours

Colour of product transparent
Closure material Low Density Polyethylene (LD-PE)
Colour of cap yellow
Tube material Polypropylene (PP)

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest bag
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 96000
Depth of case 348 mm
Width of case 278 mm
Height of case 140 mm
Case volume 0.0135 cbm
Weight of product 0.0016 kg
Weight of case 1.86 kg
EAN of inner box 4038917129596
EAN case 4038917022309

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT