



Product description

Order number	41.1506.005
Product description	Micro sample tube Citrate 9NC 0.106 mol/l 3.2%, venous blood collection, preparation: Citrate 9NC, nominal volume: 1.3 ml, (LxØ) with cap: 44 x 10.8 mm, push cap, cap: green, colour code EU, with plastic label, label/print: green, skirted conical base, 100 piece(s)/polystyrene container

Product characteristics

Type of collection	venous
Type of preparation	Citrate 9NC
Preparation concentration	0.106 mol/l
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	green
Graduation	yes
Cap	flat cap
Closure type	push cap
Application	venous blood collection
Base shape	skirted conical base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**Micro sample tube Citrate 9NC 0.106 mol/l 3.2%,
1.3 ml, push cap, EU**

Page 2

Volume of work	1,3 ml
Recommended centrifugation	1500 x g
Centrifugation time	10 Min

Size

Sample volume	1.3 ml
Diameter	10.8 mm
Length of product	44 mm
Length including cap	44 mm
Length excluding cap	40 mm

Material & colours

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	green
Tube material	Polypropylene (PP)

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	polystyrene container
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	44000
Depth of case	568 mm
Width of case	348 mm
Height of case	155 mm
Case volume	0.0306 cbm
Weight of product	0.0017 kg
Weight of case	2.47 kg
EAN of inner box	4038917129626
EAN case	4038917022316

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-23

This document was prepared by EDP support and is valid without signature.



SARSTEDT