



## Product description

|                     |                                                                                                                                                                                                                                                                            |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 41.3395.005                                                                                                                                                                                                                                                                |
| Product description | Micro sample tube EDTA K3E, venous blood collection, preparation: K3 EDTA, nominal volume: 1.3 ml, (LxØ) with cap: 47 x 10.8 mm, screw cap, cap: red, colour code EU, with paper label, label/print: red, with pierceable membrane, skirted conical base, 100 piece(s)/bag |

## Product characteristics

|                           |                         |
|---------------------------|-------------------------|
| Type of collection        | venous                  |
| Type of preparation       | K3 EDTA                 |
| Preparation concentration | 1.6 mg/ml blood         |
| Colour code               | EU                      |
| Label/ Print              | with paper label        |
| Colour of print/label     | red                     |
| Graduation                | yes                     |
| Cap                       | membrane cap            |
| Closure type              | screw cap               |
| Application               | venous blood collection |
| Base shape                | skirted conical base    |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Micro sample tube EDTA K3E, 1.3 ml, screw cap,  
EU**

**Page 2**

Volume of work 1,3 ml

**Size**

Sample volume 1.3 ml  
Diameter 10.8 mm  
Length of product 47 mm  
Length including cap 47 mm  
Length excluding cap 45 mm

**Material & colours**

Colour of product transparent  
Closure material Polypropylene (PP)  
Colour of cap red  
Tube material Polypropylene (PP)

**Purity & certification**

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

**Packaging**

Minimum order qty. 1000  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 100  
Piece (s) / outer case 1000  
Piece(s) / pallet 96000  
Depth of case 348 mm  
Width of case 278 mm  
Height of case 140 mm  
Case volume 0.0135 cbm  
Weight of product 0.0017 kg  
Weight of case 1.97 kg  
EAN of inner box 4038917265416  
EAN case 4038917265409

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**