



Product Specification
Sample tube, 10 ml, cap natural, (LxØ): 98.25 x 16.15 mm, with paper label

Page 1



Product description

| | |
|---------------------|--|
| Order number | 46.9924.106 |
| Product description | Sample tube, nominal volume: 10 ml, PP, cap natural, screw cap, (LxØ): 98.25 x 16.15 mm, flat base, with paper label, 100 piece(s)/bag |

Product characteristics

| | |
|----------------------------|--|
| Label/ Print | with paper label |
| Colour of print/label | black |
| Graduation | no |
| Cap | standard screw cap |
| Closure type | screw cap |
| Application | collection and storage of diagnostic specimens |
| Base shape | flat base |
| Volume of work | 10 ml |
| Centrifugation max (RCF) | 5000 x g |
| Centrifugation temperature | 21 |
| Centrifugation time | 15 Min |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Sample tube, 10 ml, cap natural, (LxØ): 98.25 x 16.15 mm, with paper label

Page 2

Size

| | |
|-------------------|----------|
| Sample volume | 10 ml |
| Diameter | 16.15 mm |
| Length of product | 98.25 mm |

Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | translucent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | natural |
| Tube material | Polypropylene (PP) |

Purity & certification

| | |
|---------------------------|-----------------------|
| Satisfies the requirement | sealed at 18"hg/10min |
| Latex-free | yes |
| Pyrogen-free | no |
| Safety data sheet | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 48000 |
| Depth of case | 604 mm |
| Width of case | 264 mm |
| Height of case | 373 mm |
| Case volume | 0.06 cbm |
| Weight of case | 7.6 kg |
| EAN of inner box | 4038917154895 |
| EAN case | 4038917036467 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

