



**Product Specification**  
**S-Monovette® Lithium heparin LH, 7.5 ml, cap green, (LxØ): 92 x 15 mm, with paper label**

**Page 1**



**Product description**

Order number	01.1604.100
Product description	S-Monovette® Lithium heparin LH, preparation: Lithium heparin, 7.5 ml, membrane screw cap, cap green, colour code ISO, (LxØ) without cap: 92 x 15 mm, with paper label, label/print: green, 50 piece(s)/case, sterile

**Product characteristics**

Type of collection	venous
Type of preparation	Lithium heparin
Preparation concentration	16 IU/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.







**Product Specification**  
**S-Monovette® Lithium heparin LH, 7.5 ml, cap green, (LxØ): 92 x 15 mm, with paper label**

Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	129 mm
Case volume	0.0227 cbm
Weight of product	0.0084 kg
Weight of case	5.05 kg
EAN of inner box	4038917131919
EAN case	4038917024037

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

