



Product Specification
S-Monovette® EDTA K3E, 7.5 ml, cap red, (LxØ):
92 x 15 mm, with plastic label

Page 1



Product description

Order number	01.1605
Product description	S-Monovette® EDTA K3E, K3E, preparation: K3 EDTA, 7.5 ml, membrane screw cap, cap red, colour code EU, (LxØ) without cap: 92 x 15 mm, with plastic label, label/print: red, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	red
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Liquid-tight	yes
Connection	membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® EDTA K3E, 7.5 ml, cap red, (LxØ):
92 x 15 mm, with plastic label

Page 2

Size

Sample volume	7.5 ml
Diameter	15 mm
Length including cap	108 mm
Length excluding cap	92 mm
Length including cap and plunger	129 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	20000
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	32000
Depth of box	189 mm
Width of box	134 mm
Height of box	82 mm
Depth of case	400 mm
Width of case	144 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
S-Monovette® EDTA K3E, 7.5 ml, cap red, (LxØ):
92 x 15 mm, with plastic label

Height of case	440 mm
Case volume	0.0253 cbm
Weight of product	0.0076 kg
Weight of case	4.65 kg
EAN of inner box	4038917557849
EAN case	4038917557832

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

