




## Product Specification

**S-Monovette® Lithium heparin gel LH, 4.7 ml, cap green, (LxØ): 75 x 15 mm, with paper label**

**Page 1**



### Product description

Order number	03.1631.100
Product description	S-Monovette® Lithium heparin gel LH, preparation: Lithium heparin gel, 4.7 ml, membrane screw cap, cap green, colour code ISO, (LxØ) without cap: 75 x 15 mm, with paper label, label/print: green, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin gel
Preparation concentration	25 IU/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	15 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Lithium heparin gel LH, 4.7 ml, cap green, (LxØ): 75 x 15 mm, with paper label**

**Page 2**

Liquid-tight Connection	yes membrane screw cap
-------------------------	---------------------------

### Size

Sample volume	4.7 ml
Diameter	15 mm
Length including cap	91 mm
Length excluding cap	75 mm
Length including cap and plunger	112 mm

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	green
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	26000
Depth of box	189 mm
Width of box	117 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® Lithium heparin gel LH, 4.7 ml, cap green, (LxØ): 75 x 15 mm, with paper label**

Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	129 mm
Case volume	0.0227 cbm
Weight of product	0.0075 kg
Weight of case	4.54 kg
EAN of inner box	4038917132213
EAN case	4038917024327

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

