




## Product Specification

**S-Monovette® Serum Gel CAT LightPROTECT, 4 ml, cap brown, (LxØ): 75 x 13 mm, with plastic label**

**Page 1**



### Product description

Order number	04.1925.003
Product description	S-Monovette® Serum Gel CAT LightPROTECT, preparation: Clotting activator/gel, 4 ml, membrane screw cap, cap brown, colour code EU/ISO, (LxØ) without cap: 75 x 13 mm, with plastic label, label/print: transparent/brown, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Clotting activator/gel
Colour code	EU/ISO
Label/ Print	with plastic label
Colour of print/label	transparent/brown
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**S-Monovette® Serum Gel CAT LightPROTECT, 4 ml, cap brown, (LxØ): 75 x 13 mm, with plastic label**

**Page 2**

Centrifugation time	10 Min
Liquid-tight	yes
Connection	membrane screw cap

### Size

Sample volume	4 ml
Diameter	13 mm
Length including cap	91 mm
Length excluding cap	75 mm
Length including cap and plunger	112 mm

### Material & colours

Colour of product	white
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	brown
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	50
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	39000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® Serum Gel CAT LightPROTECT, 4 ml, cap brown, (LxØ): 75 x 13 mm, with plastic label**

**Page 3**

Depth of box	171 mm
Width of box	116 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0054 kg
Weight of case	3.47 kg
EAN of inner box	4038917528894
EAN case	4038917528887

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**