




Product Specification
S-Monovette® EDTA K3E, 4.9 ml, cap red, (LxØ):
90 x 13 mm, with paper label

Page 1



Product description

| | |
|---------------------|---|
| Order number | 04.1931.010 |
| Product description | S-Monovette® EDTA K3E, preparation: K3 EDTA, 4.9 ml, membrane screw cap, cap red, colour code EU, (LxØ) without cap: 90 x 13 mm, with paper label "Blutbank", label/print: red, 50 piece(s)/case, sterile |

Product characteristics

| | |
|---------------------------|--------------------|
| Type of collection | venous |
| Type of preparation | K3 EDTA |
| Preparation concentration | 1.6 mg/ml blood |
| Colour code | EU |
| Label/ Print | with paper label |
| Colour of print/label | red |
| Graduation | fill mark |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Liquid-tight | yes |
| Connection | membrane screw cap |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® EDTA K3E, 4.9 ml, cap red, (LxØ):
90 x 13 mm, with paper label

Page 2

Size

| | |
|----------------------------------|--------|
| Sample volume | 4.9 ml |
| Diameter | 13 mm |
| Length including cap | 106 mm |
| Length excluding cap | 90 mm |
| Length including cap and plunger | 127 mm |

Material & colours

| | |
|-------------------------|-----------------------------------|
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | red |
| Membrane material | rubber |
| Tube material | Polypropylene (PP) |
| Piston material | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS) |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | IATA, ADR |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|--------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 36000 |
| Depth of box | 171 mm |
| Width of box | 130 mm |
| Height of box | 73 mm |
| Depth of case | 388 mm |
| Width of case | 358 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® EDTA K3E, 4.9 ml, cap red, (LxØ):
90 x 13 mm, with paper label

| | |
|-------------------|---------------|
| Height of case | 150 mm |
| Case volume | 0.0208 cbm |
| Weight of product | 0.0049 kg |
| Weight of case | 3.19 kg |
| EAN of inner box | 4038917134507 |
| EAN case | 4038917077545 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

