



Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 3 ml,
cap blue, (LxØ): 66 x 11 mm, with paper label**

Page 1



Product description

| | |
|---------------------|---|
| Order number | 05.1165.100 |
| Product description | S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, preparation: Citrate 9NC, 3 ml, membrane screw cap, cap blue, colour code ISO, (LxØ) without cap: 66 x 11 mm, with paper label, label/print: blue, 50 piece(s)/case, sterile |

Product characteristics

| | |
|----------------------------|--------------------|
| Type of collection | venous |
| Type of preparation | Citrate 9NC |
| Preparation concentration | 0.106 mol/l |
| Colour code | ISO |
| Label/ Print | with paper label |
| Colour of print/label | blue |
| Graduation | fill mark |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Recommended centrifugation | 1800 x g |
| Centrifugation temperature | 20 |
| Centrifugation time | 10 Min |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 3 ml,
cap blue, (LxØ): 66 x 11 mm, with paper label**

Page 2

| | |
|-------------------------|---------------------------|
| Liquid-tight Connection | yes membrane screw cap |
|-------------------------|---------------------------|

Size

| | |
|----------------------------------|--------|
| Sample volume | 3 ml |
| Diameter | 11 mm |
| Length including cap | 82 mm |
| Length excluding cap | 66 mm |
| Length including cap and plunger | 103 mm |

Material & colours

| | |
|-------------------------|-----------------------------------|
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | blue |
| Membrane material | rubber |
| Tube material | Polypropylene (PP) |
| Piston material | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS) |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | IATA, ADR |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|--------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 56000 |
| Depth of box | 145 mm |
| Width of box | 108 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 3 ml,
cap blue, (LxØ): 66 x 11 mm, with paper label**

| | |
|-------------------|---------------|
| Height of box | 63 mm |
| Depth of case | 333 mm |
| Width of case | 303 mm |
| Height of case | 125 mm |
| Case volume | 0.0126 cbm |
| Weight of product | 0.0041 kg |
| Weight of case | 2.61 kg |
| EAN of inner box | 4038917126113 |
| EAN case | 4038917020442 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

