



Product description

| | |
|---------------------|---|
| Order number | 01.2040.001 |
| Product description | S-Monovette® cfDNA Exact, preparation: cfDNA stabiliser, 9.2 ml, membrane screw cap, cap raspberry coloured, (LxØ) without cap: 100 x 15 mm, with paper label, label/print: dark red, 20 piece(s)/case, sterile |

Product characteristics

| | |
|-----------------------|--------------------|
| Type of collection | venous |
| Type of preparation | cfDNA stabiliser |
| Label/ Print | with paper label |
| Colour of print/label | dark red |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Liquid-tight | yes |
| Connection | membrane screw cap |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

S-Monovette® cfDNA Exact, 9.2 ml, cap raspberry coloured, (LxØ): 100 x 15 mm, with paper label

Page 2

Size

| | |
|----------------------------------|--------|
| Sample volume | 9.2 ml |
| Diameter | 15 mm |
| Length including cap | 116 mm |
| Length excluding cap | 100 mm |
| Length including cap and plunger | 136 mm |

Material & colours

| | |
|-------------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | raspberry coloured |
| Membrane material | rubber |
| Tube material | Polypropylene (PP) |
| Piston material | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS) |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | IATA, ADR |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|--------|
| Minimum order qty. | 80 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 20 |
| Piece (s) / outer case | 80 |
| Piece(s) / pallet | 8640 |
| Depth of box | 184 mm |
| Width of box | 141 mm |
| Height of box | 83 mm |
| Depth of case | 388 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® cfDNA Exact, 9.2 ml, cap raspberry coloured, (LxØ): 100 x 15 mm, with paper label

| | |
|-------------------|---------------|
| Width of case | 345 mm |
| Height of case | 105 mm |
| Case volume | 0.0141 cbm |
| Weight of product | 0.0099 kg |
| Weight of case | 1.9 kg |
| EAN of inner box | 4038917576666 |
| EAN case | 4038917576659 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

