



## Product description

|                     |                                                                                                                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 02.1726.021                                                                                                                                                                                 |
| Product description | S-Monovette® neutral Z, without preparation, 9 ml,<br>membrane screw cap, cap white, (LxØ) without cap: 92 x 16<br>mm, with paper label, label/print: white, 1 piece(s)/blister,<br>sterile |

## Product characteristics

|                            |                     |
|----------------------------|---------------------|
| Type of collection         | venous              |
| Type of preparation        | without preparation |
| Label/ Print               | with paper label    |
| Colour of print/label      | white               |
| Graduation                 | fill mark           |
| Cap                        | membrane screw cap  |
| Closure type               | screw cap           |
| Recommended centrifugation | 2000 x g            |
| Centrifugation temperature | 20                  |
| Centrifugation time        | 10 Min              |
| Liquid-tight               | yes                 |
| Connection                 | membrane screw cap  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Size

|                                     |        |
|-------------------------------------|--------|
| Sample volume                       | 9 ml   |
| Diameter                            | 16 mm  |
| Length including cap                | 108 mm |
| Length excluding cap                | 92 mm  |
| Length including cap and<br>plunger | 129 mm |

## Material & colours

|                         |                                   |
|-------------------------|-----------------------------------|
| Colour of product       | transparent                       |
| Closure material        | High Density Polyethylene (HD-PE) |
| Colour of cap           | white                             |
| Membrane material       | rubber                            |
| Tube material           | Polypropylene (PP)                |
| Piston material         | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS)                  |

## Purity & certification

|                           |                                |
|---------------------------|--------------------------------|
| Satisfies the requirement | IATA, ADR                      |
| Product category          | In vitro diagnostic device, CE |
| CE certified              | CE                             |
| Purity standard           | sterile                        |
| Sterilisation             | Electron irradiation           |
| Batched                   | yes                            |

## Packaging

|                                  |         |
|----------------------------------|---------|
| Minimum order qty.               | 500     |
| Type of smallest<br>subpackaging | blister |
| Piece(s) / inner box             | 100     |
| Piece (s) / outer case           | 500     |
| Piece(s) / pallet                | 13500   |
| Depth of box                     | 351 mm  |
| Width of box                     | 151 mm  |
| Height of box                    | 171 mm  |
| Depth of case                    | 778 mm  |
| Width of case                    | 368 mm  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® neutral Z, 9 ml, cap white, (LxØ): 92  
x 16 mm, with paper label**

|                   |               |
|-------------------|---------------|
| Height of case    | 195 mm        |
| Case volume       | 0.0558 cbm    |
| Weight of product | 0.0078 kg     |
| Weight of case    | 5.5 kg        |
| EAN of inner box  | 4038917133234 |
| EAN case          | 4038917024761 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

