



Product Specification
Monovette® Luer neutral Z, 9 ml, cap natural,
(LxØ): 92 x 16 mm, with paper label

Page 1



Product description

| | |
|---------------------|--|
| Order number | 02.226.160 |
| Product description | Monovette® Luer neutral Z, without preparation, nominal volume: 9 ml, (LxØ): 92 x 16 mm, Luer screw cap, cap natural, with paper label, label/print: black, animal venous blood collection, 100 piece(s)/bag |

Product characteristics

| | |
|----------------------------|--------------------------------|
| Type of collection | venous |
| Type of preparation | without preparation |
| Colour code | ISO |
| Label/ Print | with paper label |
| Colour of print/label | black |
| Cap | Luer screw cap |
| Closure type | screw cap |
| Application | animal venous blood collection |
| Recommended centrifugation | 2000 x g |
| Centrifugation time | 10 Min |
| Connection | Luer |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Monovette® Luer neutral Z, 9 ml, cap natural,
(LxØ): 92 x 16 mm, with paper label

Page 2

Size

| | |
|----------------------------------|--------|
| Sample volume | 9 ml |
| Diameter | 16 mm |
| Height of product | 92 mm |
| Length of product | 92 mm |
| Length including cap | 117 mm |
| Length excluding cap | 92 mm |
| Length including cap and plunger | 127 mm |

Material & colours

| | |
|-------------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | natural |
| Tube material | Polypropylene (PP) |
| Piston material | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS) |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 100 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 20000 |
| Case volume | 0.0382 cbm |
| Weight of product | 0.0101 kg |
| Weight of case | 5.3 kg |
| EAN of inner box | 4038917104326 |
| EAN case | 4038917032254 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-25

This document was prepared by EDP support and is valid without signature.



SARSTEDT