



## Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 1.8 ml, cap green, (LxØ): 75 x 13 mm, with paper label**

**Page 1**



### Product description

Order number	04.1955.001
Product description	S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, preparation: Citrate 9NC, 1.8 ml, membrane screw cap, cap green, colour code EU, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: green, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Citrate 9NC
Preparation concentration	0.106 mol/l
Colour code	EU
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	3000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 1.8 ml, cap green, (LxØ): 75 x 13 mm, with paper label**

**Page 2**

Liquid-tight Connection      yes  
membrane screw cap

### Size

Sample volume      1.8 ml  
Diameter      13 mm  
Length including cap      91 mm  
Length excluding cap      75 mm  
Length including cap and plunger      112 mm

### Material & colours

Colour of product      transparent  
Closure material      High Density Polyethylene (HD-PE)  
Colour of cap      green  
Membrane material      rubber  
Tube material      Polypropylene (PP)  
Piston material      High Density Polyethylene (HD-PE)  
Piston plunger material      Polystyrene (PS)

### Purity & certification

Satisfies the requirement      IATA, ADR  
Product category      In vitro diagnostic device, CE  
CE certified      CE  
Purity standard      sterile  
Sterilisation      Electron irradiation  
Batched      yes

### Packaging

Minimum order qty.      500  
Type of smallest subpackaging      case  
Piece(s) / inner box      50  
Piece (s) / outer case      500  
Piece(s) / pallet      39000  
Depth of box      171 mm  
Width of box      116 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 1.8 ml, cap green, (LxØ): 75 x 13 mm, with paper label**

**Page 3**

Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0049 kg
Weight of case	3.18 kg
EAN of inner box	4038917376808
EAN case	4038917376792

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**