




**Product Specification**  
**Urine measuring system, 2 l, with anti-reflux valve, exchange bag 2,000 ml, sterile**

**Page 1**



**Product description**

Order number	74.5681.200
Product description	Urine measuring system, 2 l, with anti-reflux valve, exchange bag 2,000 ml, sterile, tube length: 120 cm, 1 piece(s)/blister

**Product characteristics**

Colour of print/label	blue
Design	exchange bag 2,000 ml
Application	Urine drainage
Volume of work	500 ml
Anti-reflux valve	yes

**Size**

Sample volume	2 l
Max. volume	2 l
Diameter of tube	10.1 mm
Length of tube	120 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine measuring system, 2 l, with anti-reflux  
valve, exchange bag 2,000 ml, sterile**

**Page 2**

**Material & colours**

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

**Purity & certification**

Product category	Medical device
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

**Packaging**

Minimum order qty.	10
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	10
Piece(s) / pallet	280
Depth of case	580 mm
Width of case	290 mm
Height of case	235 mm
Case volume	0.0395 cbm
Weight of product	0.3364 kg
Weight of case	4.04 kg
EAN of inner box	6955824024425
EAN case	6955824024043

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-06

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**