



Product Specification
Urine bag, 2 l, with anti-reflux valve, with drainage valve

Page 1



Product description

Order number	74.5228.101
Product description	Urine bag, 2 l, with anti-reflux valve, with drainage valve, tube length: 90 cm, 150 piece(s)/case
PZN number	16333732

Product characteristics

Colour of print/label	blue
Design	with drainage valve
Application	Urine drainage
Volume of work	2 l
Anti-reflux valve	yes

Size

Sample volume	2 l
Max. volume	2 l
Diameter of tube	6.6 mm
Length of tube	90 cm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polyvinylchloride (PVC)
Colour of product	transparent

Purity & certification

Product category	Medical device
CE certified	CE supplier
Batched	yes

Packaging

Minimum order qty.	150
Type of smallest subpackaging	case
Piece(s) / inner box	150
Piece (s) / outer case	150
Piece(s) / pallet	6600
Depth of case	580 mm
Width of case	390 mm
Height of case	150 mm
Case volume	0.0339 cbm
Weight of product	0.0509 kg
Weight of case	7.94 kg
EAN of inner box	6971243660990
EAN case	6971243660648

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

