



Product description

| | |
|---------------------|-----------------------------------|
| Order number | 74.5767.003 |
| Product description | UroStretch, size: L, material: PL |
| Aid number | 15.25.11.1034 |
| PZN number | 8406960 |

Material & colours

| | |
|-------------------|----------------|
| Product material | polyester (PL) |
| Colour of product | white |

Purity & certification

| | |
|------------------|----------------|
| Product category | Medical device |
| CE certified | CE supplier |
| Batched | yes |

Packaging

| | |
|-------------------------------|------|
| Minimum order qty. | 10 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 10 |
| Piece (s) / outer case | 10 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
UroStretch, size: L

Page 2

| | |
|-------------------|---------------|
| Piece(s) / pallet | 2340 |
| Depth of case | 283 mm |
| Width of case | 173 mm |
| Height of case | 135 mm |
| Case volume | 0.0066 cbm |
| Weight of product | 0.043 kg |
| Weight of case | 0.55 kg |
| EAN of inner box | 5703936100121 |
| EAN case | 5703936100169 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT