



Product description

| | |
|---------------------|--|
| Order number | 79.8710.576 |
| Product description | Ink ribbon, Consumables for protocol printer SP542MD, violet |

Product characteristics

| | |
|-------------|---------|
| Application | Warming |
|-------------|---------|

Size

| | |
|-------------------|-------|
| Width of product | 7 cm |
| Length of product | 12 cm |

Material & colours

| | |
|-------------------|--------|
| Colour of product | violet |
|-------------------|--------|

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 1 |
| Piece(s) / pallet | 32 |
| Depth of case | 440 mm |
| Width of case | 390 mm |
| Height of case | 220 mm |
| Case volume | 0.0378 cbm |
| Weight of product | 3.7 kg |
| Weight of case | 4.5 kg |
| EAN of inner box | 4038917975902 |
| EAN case | 4038917975919 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

