



Product description

Order number	79.8710.577
Product description	Ink ribbon, Consumables for protocol printer SP742MD, black

Product characteristics

Application	Warming
Suitable for	Protocol printer

Size

Width of product	7 cm
Height of product	3 cm
Length of product	12 cm

Material & colours

Colour of product	black
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	32
Depth of case	440 mm
Width of case	390 mm
Height of case	220 mm
Case volume	0.0378 cbm
Weight of product	3.7 kg
Weight of case	4.5 kg
EAN of inner box	4038917975940
EAN case	4038917975957

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

