



## **Product description**

Order number	97.8710.541
Product description	Service module SAHARA, for adjusting and calibrating the SAHARA-III, TRANSTHERM-III and SAHARA-TSC

## **Product characteristics**

Application	Warming
-------------	---------

### **Size**

Width of product	28 cm
Height of product	7.5 cm
Length of product	33 cm

## **Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Case volume	0.0314 cbm
Weight of product	1.85 kg
Weight of case	2.18 kg
EAN of inner box	4028339753226
EAN case	4028339753219

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

