



Product Specification
**MEDISPONGE, impregnated with 4%
chlorhexidine gluconate solution**

Page 1



Product description

Order number	54.9912.100
Product description	MEDISPONGE, surgical scrub sponge for pre-operative cleaning of hands and fingernails, impregnated with 4% chlorhexidine gluconate solution, material: PUR, (LxWxH): 130 x 63 x 30 mm

Product characteristics

Preparation concentration	4% chlorhexidine gluconate solution
Application	surgical scrub sponge for pre-operative cleaning of hands and fingernails

Size

Width of product	63 mm
Height of product	30 mm
Length of product	130 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**MEDISPONGE, impregnated with 4%
chlorhexidine gluconate solution**

Page 2

Material & colours

Product material	Polyurethane (PUR)
Colour of product	white

Purity & certification

Latex-free	yes
Pyrogen-free	no
Safety data sheet	yes

Packaging

Minimum order qty.	450
Type of smallest subpackaging	case
Piece(s) / inner box	75
Piece (s) / outer case	450
Piece(s) / pallet	7200
Depth of case	528 mm
Width of case	465 mm
Height of case	526 mm
Case volume	0.1283 cbm
Weight of case	7.9 kg
EAN of inner box	4038917150491
EAN case	4038917036856

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-06

This document was prepared by EDP support and is valid without signature.



SARSTEDT