



### **Product description**

Order number	74.4701
Product description	Navel clamp, 250 piece(s)/carton
PZN number	6646559

### **Material & colours**

Product material	Polyoxymethylene (POM)
Colour of product	white

### **Purity & certification**

Product category	Medical device
CE certified	CE 0124 (DEKRA)
Purity standard	sterile
Sterilisation	Gas sterilisation
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Navel clamp, 250 piece(s)/carton**

**Page 2**

**Packaging**

Minimum order qty.	250
Type of smallest subpackaging	blister
Piece(s) / inner box	10
Piece (s) / outer case	250
Piece(s) / pallet	22500
Depth of case	388 mm
Width of case	253 mm
Height of case	165 mm
Case volume	0.0162 cbm
Weight of product	0.0064 kg
Weight of case	1.86 kg
EAN of inner box	4038917088534
EAN case	4038917006507

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-05

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**