



**Product Specification**  
**Centrifuge, supply voltage: 115 V**

**Page 1**

**Product description**

Order number 90.184.701  
Product description Centrifuge, supply voltage: 115 V

**Product characteristics**

Supply voltage 115 V

**Purity & certification**

Batched no

**Packaging**

Piece(s) / inner box 1  
Piece (s) / outer case 1  
Piece(s) / pallet 1  
Weight of product 0.0005 kg  
EAN of inner box 4038917102513  
EAN case 4038917072779

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-06

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**