



Product Specification
miniPERM® , Screw cap, for miniPERM®
production module

Page 1



Product description

| | |
|---------------------|--|
| Order number | 94.6001.036 |
| Product description | miniPERM® , Screw cap, for miniPERM® production module, not autoclavable, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister |

Product characteristics

| | |
|--------------|-----------------------------|
| Suitable for | miniPERM® production module |
|--------------|-----------------------------|

Material & colours

| | |
|------------------|--------------------|
| Product material | Polypropylene (PP) |
|------------------|--------------------|

Purity & certification

| | |
|------------------|--|
| Product category | no medical device IVD |
| Purity standard | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation | Gamma irradiation |
| Pyrogen-free | yes |
| Non-cytotoxic | yes |
| Batched | yes |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
miniPERM® , Screw cap, for miniPERM®
production module

Page 2

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 6 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 6 |
| Piece (s) / outer case | 90 |
| Piece(s) / pallet | 2160 |
| Depth of case | 283 mm |
| Width of case | 268 mm |
| Height of case | 210 mm |
| Case volume | 0.0159 cbm |
| Weight of product | 0.0013 kg |
| Weight of case | 0.91 kg |
| EAN of inner box | 4038917254724 |
| EAN case | 4038917254717 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT