




Product description

Order number	94.6001.059
Product description	miniPERM® , classic bioreactor, 2-compartment bioreactor, for protein, virus and biomass production, for suspension cells, volume of product module: 35 ml, not autoclavable, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

Product characteristics

Closure type	screw cap
Design	for suspension cells
Application	for protein, virus and biomass production
Cell type	suspension cells
Volume of work	350 ml - 450 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**miniPERM® , classic bioreactor, 2-compartment
bioreactor, for suspension cells**

Page 2

Size

Sample volume	400 ml
Production module volume	35 ml
Diameter	102 mm
Height of product	157 mm

Material & colours

Product material	Polycarbonate (PC)
Closure material	High Density Polyethylene (HD-PE)
Dialysis membrane material	Cellulose acetate (CA)

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

Packaging

Minimum order qty.	2
Type of smallest subpackaging	bag
Piece(s) / inner box	2
Piece (s) / outer case	12
Piece(s) / pallet	240
Depth of case	498 mm
Width of case	364 mm
Height of case	310 mm
Case volume	0.0562 cbm
Weight of product	0.307 kg
Weight of case	5.42 kg
EAN of inner box	4038917254847
EAN case	4038917254830

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

