




Product Specification
**Universal turning device, for 4 miniPERM®
bioreactors, 115/230 V**

Page 1



Product description

Order number	94.6001.061
Product description	Universal turning device, for 4 miniPERM® bioreactors, not autoclavable, 115/230 V, 1 piece(s)/case

Product characteristics

Suitable for	4 miniPERM® bioreactors
--------------	-------------------------

Size

Width of product	26 cm
Length of product	46 cm

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Universal turning device, for 4 miniPERM®
bioreactors, 115/230 V

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	30
Depth of case	398 mm
Width of case	358 mm
Height of case	325 mm
Case volume	0.0463 cbm
Weight of product	5.861 kg
Weight of case	6.3 kg
EAN of inner box	4038917254878
EAN case	4038917254861

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

