



**Product description**

Order number	94.6077.037
Product description	miniPERM® , Screw cap, for miniPERM® nutrient module, not autoclavable, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

**Product characteristics**

Suitable for	miniPERM® nutrient module
--------------	---------------------------

**Material & colours**

Product material	High Density Polyethylene (HD-PE)
------------------	-----------------------------------

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	16
Type of smallest subpackaging	blister
Piece(s) / inner box	16
Piece (s) / outer case	64
Piece(s) / pallet	2304
Depth of case	388 mm
Width of case	253 mm
Height of case	225 mm
Case volume	0.0221 cbm
Weight of product	0.0081 kg
Weight of case	1.2 kg
EAN of inner box	4038917255776
EAN case	4038917255769

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

