



## Product description

Order number	83.3941.102
Product description	Filtropur BT 50, Bottle top filter, for screw cap 45 mm, 500 ml, membrane: PES, Ø membraner: 90 mm, pore size: 0.1 µm, for removing mycoplasmas, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

## Product characteristics

Graduation	yes
Example of use	for removing mycoplasmas
Suitable for	screw cap 45 mm
Filter	yes
Surface area of filter	63.6 cm <sup>2</sup>
Pore size	0,1
Connection	Step taper

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Sample volume	500 ml
Filtration volume	500 ml
Height of product	10.61 cm

### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Membrane material	Polyethersulfone (PES)

### **Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

### **Packaging**

Minimum order qty.	24
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	24
Piece(s) / pallet	576
Depth of case	521 mm
Width of case	267 mm
Height of case	381 mm
Case volume	0.0529 cbm
Weight of product	0.082 kg
Weight of case	3.49 kg
EAN of inner box	4038917412292
EAN case	4038917412285

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

