



**Product description**

Order number	83.3942.100
Product description	Filtropur BT 100, Bottle top filter, for screw cap 45 mm, 1,000 ml, membrane: PES, Ø membraner: 90 mm, pore size: 0.45 µm, for clear filtration, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

**Product characteristics**

Graduation	yes
Example of use	for clear filtration
Suitable for	screw cap 45 mm
Filter	yes
Surface area of filter	63.6 cm²
Pore size	0,45
Connection	Step taper

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Size**

Sample volume	1.000.00 ml
Filtration volume	1000 ml
Height of product	15.27 cm

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent
Membrane material	Polyethersulfone (PES)

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Non-cytotoxic	yes
Batched	yes

**Packaging**

Minimum order qty.	24
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	24
Piece(s) / pallet	288
Depth of case	572 mm
Width of case	343 mm
Height of case	439 mm
Case volume	0.0861 cbm
Weight of product	0.1305 kg
Weight of case	5.08 kg
EAN of inner box	4038917412322
EAN case	4038917412315

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

